ABMA

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ABMA

Leakage Current Measuring Adapter

3-447-029-15
1/9.18

1 Applications

The ABMA leakage current measuring adapter (article number Z237A) serves as an attachment for TRMS multimeters for the measurement of touch voltage per VDE 100-710 (replaces DIN VDE 0107, passage 10), and for the measurement of continuously flowing leakage current and patient auxiliary current per DIN EN 62353:2014 and IEC 60601-1:2005/A1:2012 (Medical electrical equipment – Part 1: General requirements for basic safety and essential performance).

2 Safety Precautions

The leakage current measuring adapter includes a measuring device in accordance with DIN EN 62353 (VDE 0751-1):2015-10, figure C1. It may only be used for its intended purpose!

- Observe permissible values for continuously flowing leakage current and patient auxiliary current.
- Check the adapter for correct functioning before beginning and after finishing measurement, for example by measuring ohmic input and output resistance:

Input resistance	1 kΩ 1.0% + PTC 75 Ω
Output resistance	10 kΩ, with open input: 11 kΩ at the output terminals
Measuring shunt	1 kΩ 1.0%


Figure 1: Measuring circuit

Due to its mechanical layout, this adapter is not equipped with a fuse. Overload protection is implemented by means of a PTC element which additionally reduces the flow of leakage current in the event of a measuring point with low impedance. The impedance of the device under test should be greater than 30 kΩ. When overload protection is triggered, it's reset automatically after disconnection from the input circuit. Allow the adapter to cool down for 1 to 2 minutes in the de-energized state to this end.

! Attention!

In any case, avoid applying 230 V mains power or larger voltages to the input sockets. Although the safety of the user and downstream measuring technology is assured in any case, this may nevertheless result in impairment of the adapter's accuracy or in damage to the adapter.

Meanings of Symbols on the Instrument

Warning concerning a point of danger (attention, observe documentation!)

CE conformity marking

CAT III 300 V Measuring category III device

This device may not be disposed of with the trash. Further information regarding the WEEE mark can be accessed on the Internet by entering the search term 'WEEE' at www.gossenmetrawatt.com

3 Measuring Instrument Requirements

You'll need a measuring instrument that's capable of measuring TRMS values for direct and alternating voltages of up to 50 V. The lowest measuring range should have a resolution of at least 0.1 mV, so that the smallest permissible values for leakage current and patient auxiliary current, namely 10 μA (per DIN EN 60601-1 VDE 0750-1:2013-12), can still be detected with adequate accuracy. A measuring current of 10 μA corresponds to a display value of 10 mV.

Distance between the measuring instrument's input sockets should demonstrate the standardized value of 19 mm in order to ensure that the leakage current measuring adapter can be readily plugged in.

TRMS multimeters from GOSSEN METRAWATT are suitable for use as measuring instruments.

4 Inspecting the Leakage Current Measuring Adapter

Independent of periodic calibrations, the leakage current measuring adapter should be inspected under the following circumstances:

- after long periods of storage under unfavourable conditions,
- after overstressing of the adapter or operation beyond its specifications,
- if malfunction of the adapter is suspected (for example, if no voltage or inexplicably high voltages are indicated).

You can easily perform this inspection yourself by means of a multimeter with resistance measuring range. Please proceed as follows:

- Connect the leakage current measuring adapter directly with the resistance measurement sockets of your multimeter and measure the ohmic resistance. It should amount to 11.0 kΩ ±1%, plus the measuring uncertainty of your multimeter, if applicable. In this case, no device may be connected to the input terminals of the leakage current measuring adapter.
- Unplug the leakage current measuring adapter afterwards and connect the test cables with the resistance measurement sockets of your multimeter. Plug the test probes into the input terminals of the leakage current measuring adapter. The display value should range between 1.0 kΩ and 1.15 kΩ.

If both measurements have been completed successfully, you can use the leakage current measuring adapter for your measuring tasks. If not, further use is not permissible.

5 Technical Data

Input

Input current	0 ... 50 mA AC, DC
Measuring shunt	1 kΩ 1.0%
Overload protection	PTC
Max. input current	150 mA AC, DC

Output

Output voltage	0 ... 50 V (1 V/mA) AC, DC
Output resistance	10 kΩ
Frequency response	Per EN 60601-1

Temperature


Ambient temperature	0 ... +40 °C
Storage temperature	-25 ... +70 °C

6 Calibration

The leakage current measuring adapter is calibrated at the end of the manufacturing process. Calibration is confirmed with the attached test report. The initial calibration date of a new adapter is identical with the date of manufacture which is printed on the label with the serial number on the side of the adapter.

Example F: 52418: the 1st figure of the manufacturing date indicates the weekday, the two following figures stand for the calendar week and the last two figures for the year of manufacture. In this case it is day 5 in calendar week 24 in 2018.

Example Q28 030: Q28 is a personal inspection number and the following 3-digit number is a consecutive serial number.

7 Product Support

Technical queries (use, operation).
If required please contact:

GMC-I Messtechnik GmbH
Product Support Hotline
Phone: +49-911-8602-0
Fax: +49 911 8602-709
e-mail support@gossenmetrawatt.com

8 Repair and Replacement Parts Service Calibration Center and Rental Instrument Service

If required please contact:

GMC-I Service GmbH
Service Center
Beuthener Str. 41
90471 Nürnberg, Germany
Phone: +49-911-817718-0
Fax: +49-911-817718-253
e-mail: service@gossenmetrawatt.com
www.gmci-service.com

This address is only valid in Germany. Please contact our respective representatives or subsidiaries for service in other countries.